

20190911 - MorexSample2 - Run 10  
Adrian Stec

## Rapid Sequencing (SQK-RAD004) Guide Sheet

To print this document, navigate to:

<https://gitprint.com/MorrellLAB/Nanopore/blob/master/SQK-RAD004/SQK-RAD004-protocol.md> (<https://gitprint.com/MorrellLAB/Nanopore/blob/master/SQK-RAD004/SQK-RAD004-protocol.md>)

### Input DNA Quality Control

- ☒ Fill a large foam container with ice.
- ☒ Put on a pair of lab gloves.
- ☒ Fetch the size-selected DNA from inside a white foam container inside the lab refrigerator and place the tube on ice.
- ☒ If this is the first time using this DNA sample or the sample has been stored for a long period, you may want to check the quality of the input DNA using the Nanodrop. <sup>skip</sup> ONT recommends the following quality criteria:
  - OD 260/280 of 1.8 and OD 260/230 of 2.0-2.2
  - Concentration > 53.3 ng/uL (to be able to achieve 400 ng total in 7.5 uL)

### Check the flow cell quality

- ☒ Record today's date, the sample name, the run number, and your name at the top of this sheet.
  - Follow the convention "YYYYMMDD\_SampleName\_RunNumber" for the run number and sample name.
  - Example: 20190403\_MorexSample2\_Run3
- ☐ Check that the computer has at 500 GB of free space (not counting any connected external drives). <sup>100GB</sup> 665 GB available
- ☐ Fetch the flow cell from the lab refrigerator and record the flow cell information.
  - Flow cell type (FLO-MIN106 or FLO-MIN107): FLO-MIN106 R9
  - Flow cell label (Example: 11000540): 11000540
  - Flow cell barcode (Example: 2002540064): 2002525835
- ☒ Plug the MinION into the computer.
- ☒ Remove the configuration test cell and place it in the MinION box for safekeeping.
- ☒ Insert the flow cell into the MinION.

- the dock*
 ■ ☒ Click "MinKNOW" in the upper left, then "About MinKNOW", and then record the MinKNOW version numbers.
  - MinKNOW version (Example: 19.06.7): *19.06.8*
  - Core version (Example: 3.4.5): *3.4.8*
  - Bream version (Example: 4.1.8): *4.1.9*
  - GUI version (Example: 3.4.12): *3.4.5*
  - Guppy version (Example: 3.0.4): *3.0.7*
- ☒ Select the flowcell, and change the flow cell type to match the flow cell being used.
- ☒ Click "Check flow cells" at the bottom of the screen. Then click "Start test"
- ☐ Continue to library preparation while the test is running. When the test completes, make sure that the active pore count is > 800.
  - Number of active pores: *1560*
  - Flow cell ID (Example: FAJ01802): *FAH84690*

## Library preparation

- ☒ Fetch the following reagents from the freezer and record their batch numbers.
  - FRA (brown) batch number: *SK1311003*
  - RAP (green) batch number: *SK1241003*
  - SQB (red) batch number: *SK1281002*
  - LB (pink) batch number: *SK1271003*
  - FLB (blue) (use only one vial) batch number: *SK1291004*
  - FLT (purple stiped) batch number: *SK1301003*
- ☒ Thaw the SQB (red), LB (pink), FLB (blue), and FLT (purple) while keeping the other reagents on ice.
- ☒ Briefly spin down all six reagents in the centrifuge, then immediately store on ice.
- ☒ Sanitize a pair of lab scissors using ethanol, then cut the tip off of a 10 uL pipette tip.
- ☒ Add 7.5 uL of sample DNA to a 0.2 mL thin-walled PCR tube using the cut pipette. Place the DNA sample tube back in the fridge, as it will not be used in further steps.
- ☒ With a new pipette tip, mix the vial of FRA by pipetting up and down.
- ☒ Add 2.5 uL of FRA (brown) to the PCR tube with DNA. Place the vial of FRA (brown) back on ice.
- ☒ Gently mix the PCR tube by flicking, then spin down.
- ☐ Place the PCR tube in the middle of the thermocycler plate. Avoid the wells that are

- ☒ Start the "Rapido1" cycle under the lilei username.
- ☒ Incubate the tube at 30° C for 1 minute and then at 80° C for 1 minute.
- ☒ When you hear the fan kick on, check to make sure the cycle is complete (holding at 4 C). Remove the PCR tube and place it on ice.

### Adapter attachment

- ☒ Mix the vial of RAP (green) by pipetting up and down.
- ☒ Add 1 uL of RAP (green) to the PCR tube. Place the RAP (green) vial back on ice.
- ☒ Gently mix the PCR tube by flicking, then spin down.
- ☒ Wait 13 minutes at room temperature, then store the PCR tube on ice.

### Priming the flow cell

- ☒ Vortex the SQB (red) and FLB (blue) tubes, spin down, and return them to ice.
- ☒ Open the priming port (not the SpotON sample port).
- ☒ Set a 1000 uL pipette to 200 uL and insert the tip into the priming port.
- ☒ Increase the pipette volume until you can see a small volume of buffer entering the pipette tip.
- ☒ Visually check that there is continuous buffer from the priming port across the sensor array.
- ☒ Use a pipette tip to mix the FLT (purple striped) tube by pipetting up and down.
- ☒ Add 30 uL of <sup>\* used FLT from 2 different tubes</sup> FLT (purple striped) directly to the tube of FLB (blue) and mix by pipetting up and down. Return the tube of FLT (purple striped) to ice.
- ☒ Load 800 uL of this mix into the flow cell via the priming port, avoiding the introduction of air bubbles.
- ☒ Return the rest of the FLB (blue) mixture to ice.
- ☒ Wait 5 minutes.

### Loading the flow cell

- ☒ Find a LoBind microfuge tube.
- ☒ Using a pipette tip, mix the vial of SQB (red) by pipetting up and down.
- ☒ Add 34 uL of SQB (red) to the LoBind tube, then return the SQB (red) to ice.  
*Need to mix 2 tubes of SQB to get 34 uL*

☒ Add 25.5 uL LB (pink) Mix immediately before use

- ☒ Add 4.5 uL of nuclease-free water to the LoBind tube.
- ☒ Sanitize a pair of lab scissors using ethanol, then cut the tip off of a 100 ~~u~~L pipette tip.
- ☒ Using the cut pipette tip, add all (11 uL) of the DNA library (PCR tube) to the LoBind tube.
- ☒ Open the SpotON port. You will have to break the seal. Keep the removed cover in a secure location, as you will need to replace it later.
- ☒ Load 200  $\mu$ L of the priming mix (blue) into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.
- ☒ Sanitize a pair of lab scissors using ethanol, then cut the tip off of a 100 mL pipette tip.
- ☒ Using the cut pipette tip, mix the prepared library gently by pipetting up and down just prior to loading.
- ☒ Add all (75 uL) of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.
- ☒ Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid.

### Start the MinKNOW software

- ☒ Click the "New Experiment" button at the bottom left of the GUI.
- ☒ Follow the convention "YYYYMMDD\_SampleName\_RunNumber" for both the experiment name and the sample name.
  - Example: 20190403\_MorexSample2\_Run3
- ☒ Select SQK-RAD004 as the sequencing kit.
- ☒ Disable basecalling. Leave all other settings as the defaults.
- ☒ Start the run.
- ☒ Start the program "Amphetamine" and set the duration to "Indefinitely" to keep the machine awake.
- ☒ Place all reagents back in the freezer (except the DNA sample, which should already be in the fridge).
- ☒ Stick around for 10-15 minutes to make sure that reads are being collected.
- ☒ Store this document in a safe place.